COMORE HENSIVE

真材实料·全面解析

匠心工艺设计 追求精益求精

THE PURSUIT OF CRAFTSMANSHIP DESIGN EXCELLENCE

PP STORAGE CONTACT PLATE
PP储物接盘

STORAGE PLATE COVER PLATE
储物接盘盖板

ADJUSTABLE FOOT RESTRAINT
可调节脚托

ABS LOW HOOD
ABS低罩

180-DEGREE REVERSE HANDRAIL
180度反转扶手

MANUAL RAISING LEVER
手动抬头杆

SELECTED PU LEATHER
精选PU皮革



MULTIFUNCTIONAL GYNECOLOGICAL EXAMINATION BED

多功能妇科检查床

可升降结构
ELEVATING STRUCTURE

个性化设计
PERSONALIZED DESIGN

操作方便
EASY TO OPERATE



PRODUCT PARAMETERS

产品参数



DETAILS SHOW

细节展示

1 ONE

一键升降

A KEY LIFTING

通过遥控操作，在电机作用下调节高度
适应人体各个部位的弯曲度
使颈部、背部、臀部等分别处于舒适的状态

BY REMOTE CONTROL OPERATION, UNDER THE ACTION OF THE MOTOR TO ADJUST THE HEIGHT
ADAPT TO THE BENDING OF VARIOUS PARTS OF THE HUMAN BODY
MAKE THE NECK, BACK, BUTTOCKS, ETC. IN A COMFORTABLE STATE RESPECTIVELY



2 TWO

隐藏式便携托盘

HIDDEN PORTABLE TRAY

隐藏式设计，方便好用

PP塑料材质制造，耐用好清洗

HIDDEN DESIGN, EASY TO USE

PP PLASTIC MATERIAL, DURABLE AND EASY TO CLEAN



3 THREE

机械金属框架结构

HIDDEN PORTABLE TRAY

主题框架使用全金属机械结构，质感塑料外壳
承重力强，经久耐用

THE THEME FRAME USES AN ALL-METAL MECHANICAL STRUCTURE AND TEXTURED PLASTIC SHELL
STRONG BEARING CAPACITY AND DURABILITY



4 FOUR

背部/腿部气杆调节手柄

BACK/LEG AIR ROD ADJUSTING HANDLE





5 FIVE

腿托双调节设计

LEG SUPPORT DOUBLE ADJUSTMENT DESIGN



6 SIX

欧盟CE认证电机

EUROPEAN CE CERTIFIED MOTOR



7 SEVEN

PU皮革


THE PU LEATHER



8 EIGHT

高回弹海绵


HIGH REBOUND SPONGE



9 SEVEN

可拆卸式扶手


REMOVABLE HANDRAIL



10 EIGHT

配置孔位塞

EQUIPPED WITH HOLE PLUG



UP AND DOWN

功能展示·多角度调节

双马达电动多段升降 角度任意调节 适用多种行业

DOUBLE MOTOR ELECTRIC MULTI - SECTION LIFTING ANGLE ARBITRARY ADJUSTMENT FOR A VARIETY OF INDUSTRIES



抬头0-80°（手动）
HEAD UP 0-80°（MANUAL）

脚托0-45°（手动）
FOOT SUPPORT 0-45°（MANUAL）

底座可升降(电动)
ADJUSTABLE (ELECTRIC)